



Accuracy  
Reliability  
Intelligence

*Providing Quality Holter Reports to Healthcare Professionals*



## Accuracy

State of the art software

Giving you what you need to know to make an accurate diagnosis

Concise and easy to read reports

Large on-screen ECG display to ensure artifact-free setup

On-Screen map to assist with lead placement

## Reliability

On Time...Every Time

Support and training for office staff

Professional staff trained and experienced in reading and analyzing ECGs

Physician notification when dangerous arrhythmias occur

24 Hour turnaround

## Intelligence

Our monitors can record either 24 or 48 hours on a single AAA battery

Memory can not be removed to assure that ECG information will not be accidentally lost

EASI 5 Lead configuration allows us to derive a 12 lead from any strip

Files downloaded in less than 90 seconds

[www.cardiotechms.com](http://www.cardiotechms.com)

# DigiTrak XT Recorder— smart and innovative

On-screen lead map and EASI lead placement for artifact-free hookup

Lightest recorder with large onboard display, just 2.2 oz (62 g)

Recording 24 or 48 hours on a single AAA battery

EASI derived 12-lead standard on all recorders

Extra-slim design for patient comfort

Water-resistant case and 2-year warranty

Pacemaker detection with sensitivity adjustment

Swivel belt clip and pouch for easy wearability

Non-removable memory prevents accidental data loss and enhances data security




Automatic Power On self-test with onscreen status indicator guarantees reliable recording

Simple navigation buttons with tactile response for fast setup

Large, raised Event button easy for patients to find by feel

This unit is configured for:  
 First Name: JOE  
 Last Name: SMITH  
 Patient ID: 0123456789  
 Press < and ← to accept



Streamline your workflow by pre-configuring and previewing patient demographics.

Large display for each channel ensures artifact-free setups and avoids retests.

Fast downloads typically take less than 90 seconds. Our dual dock accepts DigiTrak XT and DigiTrak Plus recorders.